

FL-ThuKS30W

Laser for Kidney Stone and Prostate Treatment



FiberLAST's originally developed FL-ThuKS30W Kidney Stone and Prostate Treatment Laser was developed with the aim of providing import substitution by creating a maintenance-free and easy-to-use laser that will replace Ho:YAG lasers used in urology applications.

The laser provides an important advantage and flexibility for medical applications as it works with single phase voltage, requires no water cooling, allows 24/7 operation and is maintenance-free.



Parametreler	Değerler
Brand/Model	FiberLAST/ FL-ThuKS30W
Laser Type	Thulium Fiber Laser
Wavelength	1940 nm \pm 20 nm
Operation Mode	Pulsed and Continuous (CW)
Average Power	\geq 30W
Repeat Frequency	10 – 2400 Hz
Pulse Width	200 μ s-50 ms, adjustable

FiberLAST Fiber Lazer Sistemleri ve Teknolojileri A.Ş.

Merkez: ODTÜ Teknokent Galyum Blok Z-3, 06800 Çankaya-Ankara/ Türkiye

Şube: Mustafa Kemal Mah. Dumlupınar Blv. 280/D-E Blok E-2, METU Teknokent, 06530 Çankaya-Ankara/Türkiye

Tel +90 312 205 59 29 **Fax** +90 312 210 15 95 **E-mail** info@fiberlast.com.tr

www.fiberlast.com.tr